



3D MULTILAYER Zirconia Block

Excellent combination of aesthetics and stability



Shandong Huge Dental Material Corporation

Tel: +86 (633) 2277285 www.hugedental.com
Add: No.68 Shanhai Road, Donggang District, Rizhao City, Shandong Province, 276800, P. R. China
Fax: +86 (633) 2277298 marketing@hugedental.com





- Aesthetic effect
- Stable and durable
- Fewer steps of processes
- Multiple shades are available
- 100% Tosoh powder raw material ensure quality
- Support various systems engraving machines

3D GRADIENT

prepare for your esthetic dentistry



Seamless Gradient, Natural Aesthetics

Seamless gradient makes the color change smoothly from teeth neck to cut end. HUGE 3D Multilayer Zirconia provides you natural aesthetics of teeth.

Translucency Gradient, Bionic Transparency

From 47.5% transparency of cut end to 43% of tooth neck, unique gradient could achieve life-like texture.

Flexural Strength Gradient, Comfort and Strength

The flexural strength shows gradient from 700Mpa to 1060Mpa in order to protect opposing dentin by lower strength meanwhile ensure the strength in occlusion.



3D Multilayer Zirconia Block

could achieve a revolution of dental material and be the only Zirconia material in your lab.

